


MediByte Patient Guide

Your doctor or dentist has sent you home with the MediByte® to screen for snoring or sleep apnea.

To use the device, put the sensors on as described in this guide, slide the ON/OFF switch to the ON position (white dot), press the Event button on the front of the MediByte® and go to bed. When you wake up, press the Event button on the front of the MediByte®, remove the sensors and slide the ON/OFF switch to the OFF position.

Warnings

- Use the MediByte® as directed by your physician.
- Keep the MediByte® and its components away from flames and flammable compounds.
- Do **NOT** immerse the MediByte® in any liquids, damage will result.
- Do **NOT** bathe, shower or swim while wearing the MediByte®.
- Do **NOT** plug any MediByte® accessories into an electrical wall outlet.



1-100 Schneider Road,
Kanata, ON Canada
K2K 1Y2
+01.613.831.6690



CEpartner4U BV
Esdoornlaan 13
3951 DB Mearn